


**ETHICON**

**MONOCRYL\***

**(Poliglecaprone 25) Monofilament Suture**



# MONOCRYL\*

## (Poliglecaprone 25) Monofilament Suture

### Monofilament

For smooth passage through tissue.

### Synthetic

For minimal tissue reaction.

### Absorbable

By hydrolysis. Absorption is essentially complete between 90 and 120 days.

### Tensile Strength

DYED		UNDYED	
Post implantation	Approximate % original strength remaining	Post implantation	Approximate % original strength remaining
7 days	60%	7 days	50%
14 days	30%	14 days	20%
All of the original tensile strength is essentially lost by 28 days post-implantation.		All of the original tensile strength is essentially lost by 21 days post-implantation.	

### Colour

Violet or undyed.

### Range










Gauge sizes 6-0 to 1 (U.S.P.). Supplied as needled sutures.

### Sterilisation









Ethylene Oxide.




These sutures comply with the requirements of U.S. Pharmacopoeia (U.S.P.) and The European Pharmacopoeia (PhEUR).

# MONOCRYL\* (Poliglecaprone 25) Monofilament Suture









Needle length and description	Length	Colour	Gauge				
<b>⊙ 1/2 circle ETHIGUARD* Blunt Point</b>			2-0 3	0 3.5	1 4	2 5	4 6
MHB-1 	31mm	70cm	violet	<b>W3284</b>			
CTB 	40mm	90cm	violet		<b>W3288</b>	<b>W3289</b>	
BT-3 	45mm	90cm	violet			<b>W3298</b>	
CTXB 	48mm	90cm	violet			<b>W3291</b>	
XLH-B 	70mm	90cm	violet			<b>W3709</b>	
<b>⊙ 3/8 circle Taper Point</b>			4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
BB 	17mm	70cm	violet	<b>W3548</b>			
TE 	31mm	70cm	violet			<b>W3416</b>	
<b>⊖ 1/2 circle Taper Point plus</b>			6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
RB-2 	13mm	45cm	violet	<b>W3224</b>			
RB-1 	17mm	70cm 70cm	violet undyed			<b>W3435</b>	<b>W3437</b> <b>T4215AH</b>

# MONOCRYL\* (Poliglecaprone 25) Monofilament Suture

Needle length and description		Length	Colour	Gauge						
⊖ 1/2 circle Taper Point plus				6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
SH-1 plus 	22mm	70cm	violet			W3660	W3661	W3662		
JB VISI-BLACK* 	26mm	70cm	violet				W3664	W3665		
JB VISI-BLACK* 	26mm	8x70cm CR	violet				Y7840G			
SH plus 	26mm	8x45cm CR	violet				Y864G			
SH plus 	26mm	70cm 70cm	violet undyed				W3447 W3446	W3448		
SH plus 	26mm	70cm	violet				W3627	W3628		
MH-1 plus 	31mm	70cm 70cm	violet undyed				W3200	W3441 W3431	W3442	W3755
MH plus 	36mm	70cm	violet					W3463	W3443	W3457





⊙ 1/2 circle Taper Point (heavier wire diameters)				6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
CT 	40mm	70cm	violet						W3758	W3759
MO-45 	45mm	90cm	violet							W3727
CTX 	48mm	70cm	violet					W3730		

# MONOCRYL\* (Poliglecaprone 25) Monofilament Suture

Needle length and description		Length	Colour	Gauge						
⊙ 5/8 circle Taper Point				6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
	26mm	70cm 70cm	undyed violet				<b>T4290AH</b> <b>W3625</b>			
	26mm	70cm	violet				<b>W3637</b>	<b>A1H</b>		
	36mm	70cm	violet					<b>Y375H</b>	<b>Y376H</b>	
	36mm	70cm	violet					<b>W3604</b>		
▼ 1/2 circle TAPERCUT*				6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
	26mm	70cm 70cm	violet undyed				<b>W3439</b>	<b>W3440</b> <b>W3432</b>		
	36mm	70cm	violet					<b>W3488</b>	<b>W3489</b>	<b>W3490</b>
▼ 1/2 circle Reverse Cutting				6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3	0 3.5	1 4
	36mm	70cm	violet					<b>W3600</b>	<b>W3601</b>	
	40mm	90cm	violet							<b>W3770</b>







# MONOCRYL<sup>\*</sup>


## (Poliglecaprone 25) Monofilament Suture




Needle length and description	Length	Colour	Gauge				
▼ 3/8 circle Reverse Cutting - MultiPass <sup>*</sup> PRIME <sup>*</sup> Needle			6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
PS-3  16mm	45cm	undyed		<b>MPY500H</b>	<b>MPY501H</b>		
PS-2  19mm	45cm 45cm 70cm	violet undyed undyed			<b>MPY498H</b> <b>MPY496H</b> <b>MPY426H</b>	<b>MPY497H</b> <b>MPY427H</b>	
PS-1  24mm	45cm 70cm	undyed undyed			<b>MPY682H</b>	<b>MPY683H</b>	<b>MPY684H</b>
PS  26mm	70cm	undyed				<b>MPY3213H</b>	

# MONOCRYL<sup>\*</sup>


## (Poliglecaprone 25) Monofilament Suture

Needle length and description	Length	Colour	Gauge				
<b>▼ 3/8 circle Reverse Cutting - PRIME* Needle</b>			6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
P-1 	11mm	45cm	undyed	<b>W3214</b>			
P-3 	13mm	45cm	undyed	<b>W3215</b>	<b>W3203</b>		
PS-3 	16mm	45cm	undyed		<b>W3204</b>	<b>W3205</b>	
PS-2 	19mm	45cm	undyed			<b>W3206</b>	<b>W3207</b>
PS-1 	24mm	70cm	undyed				<b>W3208</b>
PS 	26mm	70cm	undyed			<b>W3212</b>	<b>W3213</b>

<b>▼ Compound Curve Reverse Cutting - PRIME* Needle</b>			6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
PS-4C 	16mm	45cm	undyed	<b>NP0127</b>	<b>NP0126</b>		

<b>▼ 3/8 circle Reverse Cutting</b>			6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
FS-3 	16mm	45cm	undyed	<b>W3221</b>			
FS-2 	19mm	45cm 70cm	undyed undyed	<b>W3209</b>	<b>W3201</b>	<b>W3202</b>	
FS 	26mm	70cm	undyed			<b>W3326</b>	<b>W3327</b>

# MONOCRYL<sup>\*</sup> (Poliglecaprone 25) Monofilament Suture

▼ Straight Cutting				6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
KS 	60mm	70cm	undyed			W3649	W3650	W3651

Needle length and description	Length	Colour	Gauge				
▼ 1/4 circle SPATULATED			6-0 0.7	5-0 1	4-0 1.5	3-0 2	2-0 3
S-14 	8mm	45cm	undyed	W3552			

## Ligatures and Non-needed Sutures

### 1.5m lengths

Gauge			violet 1.5m lengths 12 sachets per box
2-0	3		W3113
0	3.5		W3114
1	4		W3115

### LIGAPAK<sup>\*</sup> Dispenser Reel (radio opaque)

Gauge			violet 2.5m per reel 10 reels per box
3-0	2		W3020